Pocket Mask™

Mascara Pocket™

Masque Pocket™

NAME:
NOMBRE:
NOM:

DOCTOR:
MÉDICO:
MÉDECIN:

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Made in Taiwan ● Hecho en Taiwan ● Fabriqué à Taiwan

Description

The Pocket Mask is a clear, latex free, silicone mask with a 22 mm ID inlet fitting and a one-way valve outlet that allows air flow to leave the mask. The Pocket Mask is available in three convenient sizes: Infant (Orange), Small (Yellow), Medium (Blue).

Cautions

In all cases, breathing must be patient initiated and unassisted.
Any devices used with the Pocket Mask must allow for free flow of air into the mask.
The Pocket Mask must never be used as or with a resuscitation device.

Purpose

The Pocket Mask can be used with most inhalation devices having a 22 mm OD fitting where it is desirable for the patient to inhale from that device but have an alternate route of exhalation.

To Use

Connect the 22 mm ID port to an appropriate device as described previously. (Refer to the manufacturer's instructions for proper operation of that device). Place the Pocket Mask gently against the face with the one-way valve on the upper side of the mask near the nose. Exhalation may cause the mask to temporarily fog. Visually inspect the one-way valve during exhalation to assure that it opens properly.

Cleaning

Oils from the skin can, over time, reduce the flexibility of the silicone. It is recommended that the mask be cleaned after each use. To clean, hand wash with a mild detergent and warm water. Rinse well and allow to air dry. Inspect the mask, including the one-way valve, after every cleaning.

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